

Protocol for single barcode insertion (Galactose induction)

Written by Xianan Liu (Sasha Levy lab)

MATERIALS:

- Yeast strain(s): e.g. SHA185
- Plasmid library in bacteria: e.g. L139
- Yeast extract + Peptone + galactose liquid media (YPG)
 - Yeast extract: 10g/L
 - Peptone: 20g/L
 - Galactose: 20g/L
- 2x YPD liquid media
 - Yeast extract: 20g/L
 - Peptone: 40g/L
 - Glucose: 40g/L
- SC-ura plates
- Salmon Sperm DNA, sheared (10 mg/mL) (ThermoFisher cat. no. AM9680)

PROCEDURE:

Plasmid barcode library preparation:

1. Inoculate the whole frozen stock ($\sim 10^9$ cells) of a library into 1L of LB liquid media.
2. Incubate at 37°C overnight.
3. Prepare several frozen stocks from the overnight culture for future library preps.
4. Use QIAGEN Maxi prep kit to extract plasmid DNA.

Yeast culture preparation:

1. Inoculate the yeast strain by picking a single colony with a sterile toothpick into 50 ml of YPD in a 250 ml flask. Incubate at 30°C overnight on a rotary shaker.
2. Prepare 200 ml (30 x) or 400 ml (60 x) of 2x YPD broth.
3. After 12 ~16 h of growth, count cells using Coulter Count or hemacytometer and transfer the volume containing 6.25×10^8 cells (30 x) or 1.25×10^9 cells (60x) into 50 ml centrifuge tubes and pellet the cells at 3000g or 5 min. Resuspend the pellets in 20 ml pre-warmed 2x YPD broth and transfer to the culture flask. Add sufficient pre-warmed 2x YPD broth (30x -130 ml; 60x 280 ml) to bring the final titer to 5×10^6 cells/ml.
4. Incubate the flask at 30°C and 200 rpm until the cells have undergone two divisions, which may take up to 4 h.
5. Denature single-stranded carrier DNA (30x --350µl, 60x --650µl) for 5 min in 95°C heat block and chill in ice/water.
6. Make up appropriate volumes of transformation mix (see below) and keep in ice/water. This can be scaled up to 100x if needed.

(For one transformation, up to 10µg plasmid library can be used)

Ingredients	30x (ml)	60x (ml)
-------------	----------	----------

PEG 50% w/v	7.2	14.4
LiAc 1.0 M	1.08	2.16
ssDNA (10 mg/mL)	0.3	0.6
Plasmid library + water	2.22	4.44
Total volume	10.8	21.6

7. Harvest yeast cells by centrifugation at 3000g for 5 min. Use a number of 50 ml disposable conical tubes. Resuspend in $\frac{1}{2}$ culture volume of sterile water. Centrifuge at 3000g for 5 min and resuspend the pellet again in $\frac{1}{2}$ volume of sterile water and combine the pellets into one 50 ml conical tube. Centrifuge at 3000g for 5 min and discard the supernatant.
8. Pipette the transformation mix from step 6 onto the cell pellet and suspend the cells by vortex mixing vigorously for 1 min to make sure there are no big chunks of the pellet left.
9. Incubate the cell suspension at 42°C for 40~60 min depending on the yeast strain. (We usually use 50 min). Mix the contents of the tube by inversion at 5 min intervals to ensure a constant temperature throughout the mix.
10. Prepare YPG broth (30x -10 ml; 60x -20ml) by warming it to 30°C.
11. Centrifuge cell suspension at 3000g for 5 min. Pour off the transformation mix supernatant and remove the remaining liquid with a micropipettor.
12. Resuspend the cells gently in pre-warmed YPG broth. (30x -10 ml; 60x -20 ml), which will bring the titer close to the saturation density (assuming 2.5×10^8 is the saturation density).
13. Incubate the cell suspension at 30°C for ~16 hours for galactose-induced Cre recombination.
14. Determine the titer of the yeast culture using Coulter Count or hemacytometer. The density ($X \times 10^8$ /ml) is obtained. This number is used to estimate the number of generations that cells have undergone during galactose induction.
15. Centrifuge at 3000g for 5 min and pour off the YPG broth.
16. Resuspend cells in SC-ura broth or sterile water (30x -10ml; 60x -20 ml).
17. Perform a serial dilution (1:10, 1:100, 1:1000) to estimate the number of colonies. Plate the rest of the cell suspension (200 μ l for each plate) on SC-ura plates (30x ~ 50 plates; 60x ~100 plates).
18. Incubate the plates at 30°C for 2~3 days and count the number of colonies on serial dilution plates to estimate the total number of transformants **Y**.
19. To estimate the complexity of the library, divide the number of colonies (**Y**) by the number of cell divisions during galactose induction ($X/2.5$) is the estimated size of the library, which is $Y/(X/2.5)$ (assuming 2.5×10^8 is the saturation density).